

Laufzeitabelle GAMMAT® 21M

Desilliertes Wasser (ml)
als Wärmebad im 800 ml Glas

Goldanteil (g)

ECOLYT SG100 (ml)

ACTIVATOR SG100 (ml)

Laufzeit (h)

350 ml Becherglas mit Glas- Verdrängungskörper 350 ml Galvanokopf Rührer: Ø 6 x 35	230	1,20	80	1,20	7
		1,35	90	1,40	7
		1,50	100	1,60	7
		1,65	110	1,70	7
		1,80	120	1,90	7
		1,95	130	2,00	7
		2,10	140	2,20	7

350 ml Becherglas ohne Glas- Verdrängungskörper 350 ml Galvanokopf Rührer: Ø 7 x 40	250	2,25	150	2,40	7
		2,40	160	2,50	7
		2,55	170	2,70	7
		2,70	180	2,80	7
		2,85	190	3,00	7
275	3,00	200	3,20	7	
	3,15	210	3,40	7	
	3,30	220	3,50	7	
	3,45	230	3,60	7	
	3,60	240	3,80	7	
300	3,75	250	4,00	7	
	3,90	260	4,10	7	
	4,05	270	4,30	7	
	4,20	280	4,40	7	
	4,35	290	4,60	7	
325	4,50	300	4,70	7	
	4,65	310	4,80	7	
	4,80	320	5,00	7	
	4,95	330	5,20	7	
	5,10	340	5,40	7	
5,25	350	5,50	7		

800 ml Becherglas 800 ml Galvanokopf Rührer: Δ 15 x 42	6,00	400	6,00	8
	6,15	410	6,20	8
	6,30	420	6,40	8
	6,45	430	6,60	8
	6,60	440	6,80	8
	6,75	450	7,00	8
	6,90	460	7,20	8
	7,05	470	7,40	8
	7,20	480	7,60	8
	7,35	490	7,80	8
	7,50	500	8,00	8
	7,65	510	8,20	8
	7,80	520	8,40	8
	7,95	530	8,60	8
	8,10	540	8,80	8
	8,25	550	9,00	8
	8,40	560	9,20	8
8,55	570	9,40	8	
8,70	580	9,60	8	
8,85	590	9,80	8	
9,00	600	10,00	8	

Gramm Technik
Dentale Galvanotechnik



Gramm GmbH
Postfach 1124
D- 75229 Tiefenbronn · Mühlhausen
Telefon 072 34/9519-0
Telefax 072 34/9519-40
e-mail gramm.dental@t-online.de
<http://www.galvanoforming.de>